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

4) After confirming that there are no bubbles in arterial and venous headers, Stop the blood pump, turn the dialyser 180° to allow removal of bubbles from dialysate. Connect the dialysate lines to the dialyser and run dialysate at a flow of about 500ml/min under dialysate pressure of 0 mmHg through the extracorporeal circuit for about 5 minutes.

(Fig. 2)

Priming

Following the rinsing and air removal, start the blood pump, run about 500 ml of physiological saline (eventually heparinized, 2000 IU/500 ml) through the arterial line at a flow rate of about 100 ml/min; Prime the dialyser with heparinized saline and dialysate flow.

(Fig. 3)

Leakage test

-Activate the bypass valve on the machine and connect the dialyser lines to the machine.

-Prime the arterial, venous line and the dialyser completely with physiological saline by operating the blood pump, then stop the pump operation.

-Clamp the arterial line near the dialyser and the distal end of the venous line with forceps.


-Place the clamped distal end about 1m below the dialyser and remove the forceps. (This results in application of a negative pressure of about 70 mmHg to the blood compartment of the dialyser.)

Examine for any or no continuous stream of bubbles in the venous header to check for leakage from the dialyser; if observed, replace the dialyser with a new one.

Start of dialysis

1) Stop the blood flow. Prepare the blood access site and connect to the arterial line. Remove the forceps from the arterial and venous lines. While running dialysate at a flow of about 500 ml/min, operate the blood pump at a flow rate of about 100 ml/min.

2) Confirm that no air bubbles remain in venous header or venous bloodline.


(Fig. 4)

3) Fully prime the arterial and venous lines including the dialyser with blood by operating the blood pump, and then stop the pump operation. Clamp the distal end of the venous line with forceps.

4) Prepare the blood return site and connect to the venous line. After confirming that there are no bubbles in the line, remove the forceps from the line. After checking that there are no forceps on the lines and no line folding, operate the blood pump at a low flow rate. Avoid excessive pressure to the lines and the dialyser to prevent leakage from the dialyser and separation of each of the connection.

5) Increase the blood pump rate to 200 or 300 ml/min according to patient condition.

(Fig. 5)

Operations during dialysis

1) If stopping the blood pump is required during dialysis due to insufficient blood flow or other reasons, lower the Ultrafiltration rate UF to a safe minimum.

2) Set UF rate carefully to avoid excessive water removal. Reduce the blood flow rate accordingly to avoid a disequilibrium syndrome.

3) In case of blood leakage alarm, stop the dialysate supply and recover blood, then replace the dialyser with a new one.

Dialysis treatment-end and blood recovery

1) Lower the UF and TMP as low as possible. Stop the blood pump, clamp the arterial line and remove the line from the arterial blood access site; then connect the line to the physiological saline vial for blood recovery.

2) Unclamp the arterial line and run 100~200ml of physiological saline at a flow rate of about 100 ml/min to expel blood from the arterial and venous lines and the dialyser.

3) After blood recovery, discard the blood line sets and the dialyser by suitable means as medical waste.

Doc. Code	Rev. No.	Date Of Issue	Page
F72300RD72	R03	01/05/25 (16/Aug./2022)	4 of 4

contraindications

In general, dialysis disease is not sex-dependent. PS and PES haemodialyzers of meditechsys are not recommended for pediatric, newborn and toddler use. Factors such as diabetes, hypertension, old age and etc. Increase the likelihood of this disease.

Dialyser choice is according to its surface area, weight of patient and patient clinical test (toxins concentration) by physician.

There are no specific contraindications for the dialyser itself. Avoid use in patients with known hypersensitivity to particular dialyser (PS and PES). Usage of haemodialyzer is contraindicated in patients presenting with a known allergy to heparin or who have type II thrombocytopenia caused by heparin (HIT Syndrome type II).

PS (Polysulfone series) and PES (Polyethersulfone series) haemodialyzers are non-pyrogenic devices and sterilized with gamma irradiation.

Warning

- Do not use for any other purposes than haemodialysis.

- The dialyser should only be used as directed by a physician.

- If foaming, blood leakage, blood coagulation and hemolysis occur during the use of this product, take appropriate measures according to a physician's instructions.

- DO NOT REUSE THIS DEVICE.**

If reuse happened, risk clotting, severe infection and diminished performance is expected.

- If the patient exhibits any abnormal symptoms such as discomfort, pruritus, urticaria, asthmatic reaction, hypertension and/or arrhythmia during the use of this product, take appropriate measures according to a physician's instructions.

- Dialyser contact with chemicals and disinfections affect the properties of housing, polyurethane and fiber. The operating safety does not longer exist and the manufacturer's liability is excluded.

- Heparinize the extracorporeal circuit during the priming, these dosage of heparin and its mode of administration during dialysis are established by attending physician. The coagulation time should be checked regularly.

Device performance is diminished if it is used below mentioned flow rates, below mentioned pressure or in particular orientations. This instruction should be followed.

NB: Meditechsys will not be responsible for injury of patient or any person or damage to any object that is attributed to bad transport, bad storage, shipping accidents, abnormal use, user error and reasonably foreseeable misuse.

Caution

It is recommended to use of a dialysis machine equipped with a volumetric ultrafiltration unit, an

accurate UF control system and a precise weigh system for haemodialysis.

1. Caution before use

- Do not use the device if package is damaged or if protective end caps are not in place.

- Avoid air intake and contamination during rinsing and priming operations.

- Blood- side: Rinse extracorporeal circuit with sterile physiological saline (recommended 1000 ml and 500 ml of heparinized physiological saline, 2000 U/500ml) to prepare the fibers for blood flow at a flow rate of 100 ml/min.

- Perform a leak test to check the integrity of the bloodline and dialyser.

- Start dialysis mode immediately after priming operation.

2. Caution in use

- Continuously monitor the arterial and venous pressures in the bloodline and check for blood leakage during dialysis.

- Set TMP alarm (max. 500 mmHg).

3. Caution after use

- Dispose of the dialyser immediately after use by suitable means as contaminated medical waste and according to any prevailing environmental regulations.

- The weight loss should be monitored.

4. Caution for storage

- Store at 5°C~35°C avoiding exposure to direct sunlight, severe vibration and relative humidity between 25% -65%

5. Connectors

- The blood ports of haemodialyzers comply with the requirements of BS EN ISO 8637-2:2018 and BS EN ISO 8637-1:2020 . They can connect to female lure lock connectors of blood line set.


- The dialysis fluid port of haemodialyzer comply with the requirements of BS EN ISO 8637-1:2020 and BS EN ISO 8637-2:2018 . They can connect to the dialysis fluid port of standard dialysis machines.

Rinsing and air removal

- Unpack the dialyser and clip onto the holder ensuring that the venous side is directed upward and the arterial side is directed downward.

- Connect the arterial line to the physiological saline vial.

- Connect aseptically the arterial and venous line to the dialyser and remove air completely by running physiological saline (recommended volume 1000ml), at a flow rate of about 100ml/min

(Fig. 1)

Doc. Code	Rev. No.	Date Of Issue	Page
F72300RD72	R03	01/05/25 (16/Aug./2022)	3 of 4

Specifications

Low Flux Hollow Fiber Haemodialyzer														
Type	UF Coefficient (Ultrafiltration) (ml/mmHg.h)	KoA ml/min	Clearances (ml/min)										Blood Priming Volume (ml)	Surface Area (m ²)
			Urea		Creatinine		Phosphate		Vitamin B ₁₂		Inulin			
			Q _a = 200	Q _a = 300	Q _a = 200	Q _a = 300	Q _a = 200	Q _a = 300	Q _a = 200	Q _a = 300	Q _a = 200	Q _a = 300		
PSR 10 LF	7.3	698	187	238	167	201	143	161	79	84	-	-	59	1.0
PSR 13 LF	9.5	824	195	250	180	223	153	182	93	100	-	-	69	1.3
PSR 14 LF	11.9	977	196	253	184	226	157	188	96	117	-	-	74	1.4
PSR 16 LF	13.9	1240	197	274	189	243	164	196	120	135	-	-	86	1.6
PSR 18 LF	16.7	1570	198	284	189	265	174	215	139	156	-	-	105	1.8
PSR 20 LF	17.8	1610	199	285	191	268	177	220	145	161	-	-	107	2.0
PESR 10 LF	8.8	556	174	220	158	190	137	157	98	107	-	-	59	1.0
PESR 13 LF	10.9	689	185	237	172	212	152	180	117	131	-	-	71	1.3
PESR 16 LF	12.7	836	190	251	180	230	161	192	126	141	-	-	90	1.6
PESR 18 LF	17.9	1320	197	277	184	258	171	212	134	152	-	-	112	1.8
High Flux Hollow Fiber Haemodialyzer														
PSR 100 HF	34.6	930	190	258	173	217	160	199	113	127	78	84	59	1.0
PSR 130 HF	46.4	1030	195	264	180	234	174	219	130	146	93	103	69	1.3
PSR 140 HF	50.6	1190	196	272	186	247	179	230	139	159	104	113	74	1.4
PSR 160 HF	59.4	1570	197	284	193	267	188	249	153	178	117	130	86	1.6
PSR 180 HF	63.7	1660	198	286	194	276	192	259	161	190	125	141	105	1.8
PSR 200 HF	74.6	1900	199	290	196	283	195	277	168	208	131	152	107	2.0
PESR 130 HF	58.3	1140	196	270	192	255	181	236	147	168	107	116	72	1.3
PESR 160 HF	67.0	1610	197	285	193	269	189	254	160	186	120	134	89	1.6
PESR 180 HF	84.2	1770	198	288	197	279	194	262	165	194	127	144	110	1.8

Specifications and performance data at QB: 200/300 ml/min, QD: 500 ml/min, QF: 0 ml/min, T: 37°C
Performance data were measured in vitro according to standards BS EN ISO 8637-1:2020 and BS EN ISO 8637-2 :2018

Typical values obtained with an individual batch of fiber, the clinical use may determine difference in results in relation to different Ultrafiltration and measuring technique possible verification between batches of fibers.

UF measurement using Bovine/Human blood (Hct 32%, Protein 60 g/l)

The tolerance for clearance numbers is ± 10%

HF: High Flux LF: Low Flux

Read the label on the product or on the carton prior to use

	Date of Manufacture		Caution		Single use		Consult instructions for use		Keep away from sunlight
	Conform with relevant EU directives		Expiry date		Do not resterilize		Keep at 5°C ~ 35°C		Fragile, handle with care
	The fluid pathways (blood and dialysis fluid) that has been Sterilized using Gamma Irradiation		Authorized representative in the European Community		Manufacturer		Lot number		Do not use if package is damaged
	This way up		Use UF Controller		Non-Pyrogenic		Storage Relative humidity (RH) limit: 25% - 65%		Medical Device
			Catalogue Number		up to 7 boxes can be topped up				Keep dry

Direction for use

Intended Purpose

PS (Polysulfone series) and PES (Polyethersulfone series) haemodialyzer are mainly used for haemodialysis that is a form of renal replacement therapy (RRT) for End stage Renal Disease (ESRD) patient who has been inflicted with chronic kidney diseases. These devices are single use for patients with renal failure and under haemodialysis in which waste solutes from the blood are removed primarily by diffusion across a semi permeable membrane contained in a haemodialyzer in which the blood flows on one side of the membrane and dialysis fluid flowing on the other.

Doc. Code	Rev. No.	Date Of Issue	Page
F72300RD72	R03	01/05/25 (16/Aug./2022)	2 of 4

ICD GROUP

Meditechsys

Manufacture Address: Meditechsys Co.

Head Office: 5th Floor, No.44, Sa'adat Abad Street, Tehran, IRAN

P.O.Box: 1998898639

Tel: +98(21) 22074001 **Fax:** +98(21) 22073101

Factory: 3rd Tangestan, West Hafez Blvd., Eshtehard Industrial Zone, IRAN

Telefax: +98(26) 37773870

Website: www.meditechsys.com

E-mail: info@meditechsys.com

Always in All ways

Customer care service:

1st Floor, No.44, sa'adat Abad, Tehran, IRAN.

P.O.Box: 1998836785

Sales and marketing: +98(21) 22149688

Fax: +98(21) 22361067

Website: www.medwayteb.com

Sales@ medwayteb.com

Telegram.me/medway

@medwayteb

European Representative:

NephroCan S.r.l., via savoia 78 ,00198 ROMA, ITALY

Tel Office: +39 0685 2371 **Fax:** +39 06853 50187

www.NephroCan.com/EU-Rep

Note: Meditechsys is not responsible for any shipping accidents, bad storage and abnormal use.

Doc. Code	Rev. No.	Date Of Issue	Page
F72300RD72	R03	01/05/25 (16/Aug./2022)	1 of 4